

## ORDERING INFORMATION



### Conical Narrow Connection (CHC)

Diameter	Length	Ref. No.	Dimensions				
			A	B	C	D	H
Ø 3.2	8 mm	9308	Ø 3.2	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	10 mm	9300	Ø 3.2	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	11.5 mm	9301	Ø 3.2	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	13 mm	9303	Ø 3.2	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	16 mm	9306	Ø 3.2	Ø 2.9	Ø 1.5	Ø 2.5	2.1
Ø 3.5	8 mm	9328	Ø 3.5	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	10 mm	9320	Ø 3.5	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	11.5 mm	9321	Ø 3.5	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	13 mm	9323	Ø 3.5	Ø 2.9	Ø 1.5	Ø 2.5	2.1
	16 mm	9326	Ø 3.5	Ø 2.9	Ø 1.5	Ø 2.5	2.1



### INSERTION TOOLS

#### Manual



#### Motor Mount



#### Wrench



### Conical Standard Connection (CS)

Diameter	Length	Ref. No.	Dimensions				
			A	B	C	D	H
Ø 3.75	8 mm	9338	Ø 3.75	Ø 3.1	Ø 1.8	Ø 3.1	2.5
	10 mm	9330	Ø 3.75	Ø 2.9	Ø 1.5	Ø 3.1	2.5
	11.5 mm	9331	Ø 3.75	Ø 2.9	Ø 1.5	Ø 3.1	2.5
	13 mm	9333	Ø 3.75	Ø 2.9	Ø 1.5	Ø 3.1	2.5
	16 mm	9336	Ø 3.75	Ø 2.9	Ø 1.5	Ø 3.1	2.5
Ø 4.2	8 mm	9348	Ø 4.2	Ø 3.55	Ø 1.8	Ø 3.1	2.5
	10 mm	9340	Ø 4.2	Ø 3.3	Ø 1.8	Ø 3.1	2.5
	11.5 mm	9341	Ø 4.2	Ø 3.3	Ø 1.8	Ø 3.1	2.5
	13 mm	9343	Ø 4.2	Ø 3.3	Ø 1.8	Ø 3.1	2.5
	16 mm	9346	Ø 4.2	Ø 3.3	Ø 1.8	Ø 3.1	2.5
Ø 5.0	8 mm	9358	Ø 5.0	Ø 4.4	Ø 2.6	Ø 3.1	2.5
	10 mm	9350	Ø 5.0	Ø 4.1	Ø 2.3	Ø 3.1	2.5
	11.5 mm	9351	Ø 5.0	Ø 4.1	Ø 2.3	Ø 3.1	2.5
	13 mm	9353	Ø 5.0	Ø 4.1	Ø 2.3	Ø 3.1	2.5




### INSERTION TOOLS

#### Manual



#### Motor Mount






#### Wrench



# DRILLING PROTOCOL

## STEP DRILLING SEQUENCE

	Bone Type IV			Bone Type II & III			Bone Type I							
Ø 3.2	Ø 2.0			Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.0				
Ø 3.5	Ø 2.0	Ø 2.0 / Ø 2.4		Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.0	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.2			
Ø 3.75	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.2	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.2	Ø 3.2 / Ø 3.65 Cortical*		
Ø 4.2	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 2.8 / Ø 3.2	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 3.2 / Ø 3.65	Ø 2.0	Ø 2.4 / Ø 2.8		Ø 3.2 / Ø 3.65	Ø 3.65 / Ø 4.1 Cortical*	
Ø 5.0	Ø 2.0	Ø 2.4 / Ø 2.8	Ø 3.2 / Ø 3.65	Ø 2.0	Ø 2.4 / Ø 2.8	Ø 3.2 / Ø 3.65	Ø 3.65 / Ø 4.1	Ø 2.0	Ø 2.4 / Ø 2.8	Ø 3.2 / Ø 3.65	Ø 3.65 / Ø 4.1	Ø 4.1 / Ø 4.5	Ø 4.5 / Ø 4.8 Cortical*	

\* Cortical - Drill through cortical plate with the larger diameter.

## STRAIGHT DRILLING SEQUENCE

	Bone Type IV			Bone Type II & III				Bone Type I									
Ø 3.2	Ø 2.0			Ø 2.0	Ø 2.4	Ø 2.8*		Ø 2.0	Ø 2.8	Ø 3.0*							
Ø 3.5	Ø 2.0	Ø 2.4*		Ø 2.0	Ø 2.8	Ø 3.0*		Ø 2.0	Ø 2.8	Ø 3.2*							
Ø 3.75	Ø 2.0	Ø 2.4	Ø 2.8**	Ø 2.0	Ø 2.8	Ø 3.2**		Ø 2.0	Ø 2.8	Ø 3.2**	Ø 3.65 Cortical*						
Ø 4.2	Ø 2.0	Ø 2.8	Ø 3.2**	Ø 2.0	Ø 2.8	Ø 3.2	Ø 3.65**	Ø 2.0	Ø 2.8	Ø 3.2	Ø 3.65**	Ø 4.1 Cortical*					
Ø 5.0	Ø 2.0	Ø 2.8	Ø 3.2	Ø 3.65**	Ø 2.0	Ø 2.8	Ø 3.2	Ø 3.65	Ø 4.1**	Ø 2.0	Ø 2.8	Ø 3.2	Ø 3.65	Ø 4.1	Ø 4.5**	Ø 4.8 Cortical*	

\* cortical - Drill through cortical plate

\*\* 3mm shorter than implant's length. note that drill can be replaced by a corresponding step drill throughout entire implant's length.

**important:** professional considerations may be required for adaptations of the drill protocol in specific cases.