Elos Accurate® Scan Body

Instruction for use - Инструкция за
употреба – Návod k použití –
Brugsanvisning – Gebrauchsanleitung –
Kasutusjuhend – Instrucciones de uso –
Οδηγίες χρήσης – Mode d'emploi – Upute
za upotrebu – Istruzioni per l'uso –
Lietošanas norādījumi – Naudojimo
instrukcija – Használati útmutató –
Gebruiksaanwijzing – Bruksanvisning –
Instrukcja użytkowania – Instruções de uso
– Instrucțiuni de utilizare – Návod na
použitie – Navodila za uporabo –
Käyttöohjeet – Bruksanvisning – Kullanma
Talimatı

Disclaimer of liability: Elos Medtech Pinol products may only be used according to the instructions and recommendations. The user of Elos Medtech Pinols products should determine their suitability for particular patients and indications. Elos Medtech Pinol disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or damages arising out of or in connection with any errors in professional judgement or practice in the use or placement of the Elos Medtech Pinol products.

Please note that not all products are regulatory cleared/released for sale in all markets even though the language is present in this instruction.

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Elos Accurate® Scan Body

1 Intended use

The Scan Body is intended to be used together with a dental scanner to locate the exact position of an implant/analog or abutment in relation to prosthetic restorations. The scan body is intended to be mounted directly on the implant/analog or abutment using the Scan Body Driver. The Scan Body is intended to be used directly in a patient's mouth when using an Intra oral scanner, or on a model of the patients' teeth setup when using a desktop scanner.

2 Product Description

The Scan Body is manufactured from biocompatible titanium alloy (TiAl₆V₄ ELI) and polyetheretherketone (PEEK). The corresponding Scan Body Driver used to tighten the integrated screw is manufactured from biocompatible stainless steel. The product is available for a variety of implant platforms. For specific product descriptions, please refer to individual product labels.

3 Contraindications

None known

4 Warnings and precautions

- The Scan Body may only be used with the original branded implant/analog or abutment stated on the product label.
- Make sure that the seating surface, upon which the scan body is to be placed, is clean and free of residue and liquids as well as bone or soft tissue when mounting
- Since the Scan Body is small it must be handled with caution to avoid the patient from swallowing or inhaling it.
- Tighten the integrated screw must only be done slightly by hand using the Scan Body Driver (max. 5Ncm). To high tightening torque can damage the device which may cause inaccurate scan results.
- DO NOT use a Scan Body that has been damage or modified in any way. Use of a damage or modified Scan Body may cause inaccurate scan result.
- Do not try to disassemble the Scan Body.
- Do not try to rotate the Scan Body, when it is attached to the implant.
- The Scan Body may only be used with a compatible digital library file installed in the design software.
- The Scan Body must only be used and handled by dental professionals.
- Allergies to the materials or contents of the materials may very rarely occur.
- Inspect Scan Body/Scan Body Driver under normal lighting for the removal of visible soil and for damage or wear: Any Scan Body/Scan Body Driver showing visible damage or wear must be disposed of according to section 11.

5 Compatibility information

The Scan Body is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible with.

Scan Body's for intra oral use are marked IO. Scan Body's marked IO can also be used on models when using a desktop scanner.

Scan Body's marked SA are only intended for restorations seated in the conical interface of an implant

6 Cleaning and Sterilization

The Scan Body is delivered non-sterile. Prior to attachment of the Scan Body in the patient's mouth, it must be cleaned followed by sterilization. Further instruction can be found in "Cleaning and sterilization guideline" at <https://elosmedtech.com/IFU/>.

1. Clean according to the guidance in the "Cleaning and sterilization guideline".
2. Sterilize according to the guidance in the "Cleaning and sterilization guideline". Summarized below.

Note: For US: Steam Sterilization Cycle 132°C (270°F)/4min.

Procedure	Dynamic-air-removal sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134° C (273°F)
Drying time	20 min.	

Caution:

- As the Scan Body is a high precision tool, it may not be used while still warm from the autoclave, as the thermal expansion may lead to inaccurate scanning. Before using the Scan Body, it must be cooled to 18-37°C in temperature.
- Accompanying packaging material CANNOT be autoclaved unless it is explicitly stated that the material is suitable for sterilization.

7 Procedure

7.1 Scanning

1. Clean and sterilize the Scan Body according to the procedure described in section 6. Sterilization is only required for intra oral scanning.
2. Remove cover screw, healing abutment or similar from the implant or abutment upon which the Scan Body is intended to be placed.
3. Place the Scan Body onto the implant/ abutment in the patient's mouth or onto the analog when scanning a model.
4. Tighten the integrated screw slightly by hand using the Scan Body Driver.
5. Conduct the scan according to the scanner manufacturer's instructions to obtain a digital 3D model of the patient's mouth.
6. After scanning, the Scan Body is removed using the Scan Body Driver.

7.2 Design

1. Open the obtained digital 3D model in the scanners design software
2. Select the relevant digital library which depend on the prosthetic solution used in the specific case
3. Choose implant platform from the library corresponding to the scanned implant/analog or abutment.
4. Design the prosthetic restoration according to the design software's instructions. Be aware of that each specific prosthetic solution might have restricted design guidelines in accordance to maximum angulation and minimum wall thickness.

8 Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

9 Validity

Upon publication of this instruction for use, all previous versions are superseded.

10 Storage and Handling:

Recommended storage temperature range is 10-40°C (50-104°F) and the recommended relative humidity range is 0-80%

11 Disposal:

When disposed of the Elos Accurate® Scan Body may present an infection or microbial hazard and should therefore be disposed of as biological waste, according to national requirements.

12 Serious Incidents:

If any serious incidents occur in relation with this device, they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the

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country where the incident occurred. A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat. Contact Elos Medtech Pinol A/S at: complaint.empi@elosmedtech.com

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

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complaint.empi@elosmedtech.com adresinden Elos Medtech Pinol A/S ile temasa geçin